



Measured Parameters :
SpO₂ + PR

3000 h
SpO₂ Review

100
Patient ID

AV-ACE

Pulse Oximeter

AV-ACE

Pulse Oximeter

As a handheld pulse oximeter, AV-ACE is designed to serve the needs of both long-term monitoring and spot check. It's dedicated to offer reliable, accurate and sensitive measurement of blood oxygen saturation under even harsh conditions.

- LCD display with backlight control
- Dual work modes: monitoring & spot check
- Numeric display with plethysmogram display
- Trend review
- Pulse-tone modulation (Pitch Tone)
- Adjustable audio and visual alarms
- Powerful data storage capacity (up to 300 hours)
- PatientCare Viewer software for PC data management
- 4 x AA batteries for up to 48 hours of work
- Rechargeable battery for up to 30 hours of work
- Optional battery charger stand
- Automatic power-off function for power saving

Patientcare Viewer



PC Data Management



Easy Operation



Protective Cases



Unique Sensor Design



Bed Rail Mount



Battery Charger Stand



Carrying Bag



AVEUS INSTRUMENTS Sp. z o.o.

17/19/133 Domaniewska Street,
02-672 Warsaw Poland

ul. Domaniewska 17/19 lok.
133 02-672 Warsaw Poland

info@aveusinstruments.com.pl
www.aveusinstruments.com.pl