



AV-PRO XL

Patient Monitor

AV-PRO XL

Patient Monitor

Product Feature

- 15" Colour TFT Display, led backlight.
- Operating Mode Selection: Neonatal / Pediatric/Adult.
- Visual and Audio alarms, adjustable.
- Patient management, name and age and id.
- Network with Central Station Software,
- Built-in removable and rechargeable lithium battery,
- Easy operating UI with both Keypad menu.
- Powerful data storage and management, Trend graph, Trend table, etc
- Arrhythmia analysis and S-T segment analysis.
- Compatible to Neonatal, Pediatric, Adult Patients

Standard parameters :


3/5-Lead ECG, SpO2 , NIBP, RESP, 2-Temp, PR

Optional :

EtCO2 , Touchscreen , Thermal Recorder ,WLAN accessory , Nellcor-SPO2 ,2-IBP Masimo SpO2,Masimo AGM.

Suitable for

ICU/CCU, Operation Room,
Postoperative observatory,
Emergency room,
General ward, etc.



Standard Interface



Max 8 waveform display simultaneously
max 6 leads ECG waveforms display



OxyCRG interface



Big Font interface



Durable probes connector
with high quality



High quality Battery
2.5hours Battery Life



Built-in thermal printer



3 level Audible and visual alarming

Memory
96 hr at 30 second intervals for a NIBP
reading every 5 mintues
Display tabular; graphical



AV-PRO XL

Patient Monitor



AVEUS INSTRUMENTS Sp. z o.o.

17/19/133 Domaniewska Street,
02-672 Warsaw - Poland

ul. Domaniewska 17/19 lok.
133 02-672 Warszawa - Poland

info@aveusinstruments.com.pl
www.aveusinstruments.com.pl