



AV-SD3

Ultrasonic Pocket Doppler

AV-SD3

Ultrasonic Pocket Doppler

With excellent sensitivity, complete interchangeability and high durability, AV-SD3 series Ultrasonic Pocket Dopplers are ideal for routine fetal heart rate detection by the clinicians.



AV-SD3



AV-SD3-Pro



User-Friendly Design

- Compact design with easy one-hand operation
- Bright OLED screen with real-time FHR display
- Hi-Fi speaker

Interchangeable Waterproof Probe

- Plug-and-play
- 2 MHz probe with deep penetration for later pregnancy and larger patients
- 3 MHz probe detects fetal heart beat as early as 9 weeks
- Automatic power-off while putting the probe back to the rack (AV-SD3-Pro)

Work Mode

- Real-time work mode for spot check
- Record & play with audio playback and numeric display



Power Source

- Standard AA batteries for easy replacements
- Built-in Li-ion rechargeable battery ensures long working time
- Unique charger stand and wall mount design



Configuration

	AV-SD3	AV-SD3-Pro
FHR Display	√	√
Charger Stand	x	O
Record & Play	x	√
AA Batteries	O	x
Li-ion Battery	x	√

√ Standard x Not Available O Optional



AV-SD3

Ultrasonic Pocket Doppler



AVEUS INSTRUMENTS Sp. z o.o.

17/19/133 Domaniewska Street,
02-672 Warsaw - Poland

ul. Domaniewska 17/19 lok.
133 02-672 Warszawa - Poland

info@aveusinstruments.com.pl
www.aveusinstruments.com.pl