



A Slim
Monitoring
Experience

AV-Vital Plus

Vital Signs Monitor


AV-Vital Plus

Vital Signs Monitor

Smarter, Larger, yet Slimmer



360°
Alarm
Light



Cable
receiving
Design



<2.5 kg
Ultra-thin &
Light Weight



Built-in Barcode
Scanner



TEMP Options



8" Touch Screen

No-fan Design

LAN
DHCP

Wi-Fi
Built-in

USB&OTG

Barcode
Scan

Combining a stunning ultra-slim design, a large capacitive touch screen, together with a rich set of leading functions, AV-Vital Plus takes the lead and presents itself as the beginning of a new generation of vital signs monitors. Bringing in full options of SpO2, NIBP and Quick Temp/Infrared Ear Temp & differentiated work modes, AV-Vital Plus serves for both long time monitoring and spot check applications.

- **Standard Parameters**

Aveus SpO2, Aveus NIBP, PR

- **Optional Parameters**

Nellcor™ SpO2, SunTech® NIBP, Aveus Quick TEMP, Infrared Ear TEMP, Covidien Quick TEMP

Redefinition of Adaptive Work Modes

Round Mode

- Special work mode designed for daily ward rounding application. Worklist could be easily input through USB disk or downloaded from EMR system. Moreover, the color-coded patient list comes with square windows interface design, you can tell patients' status and information in a quick glance.

Monitoring Mode

- For continuous patient monitoring with real-time data, real-time alarms, and history trend.

Spot Mode

- For spot-check applications, such as out-patient, emergency, etc. Real-time vitals are acquired and displayed, no alarm will be seen or heard.



Patient management



Patient status with color code

TEMP Module Options



Infrared Ear TEMP



Aveus Quick TEMP



Covidien Quick TEMP



Rolling Stand

Non-volatile built-in Memory

Spot & Round Mode

Monitoring Mode (for each patient file)

16,000,000
Measurement

200
Alarm

1200
NIBP Measurement

120 h
Trend



AV-Vital Plus

Vital Signs Monitor



AVEUS INSTRUMENTS Sp. z o.o.

17/19/133 Domaniewska Street,
02-672 Warsaw Poland

ul. Domaniewska 17/19 lok.
133 02-672 Warsaw Poland

info@aveusinstruments.com.pl
www.aveusinstruments.com.pl