



30,000
NIBP Measurement
Review

72 h
Trend Review

800
Alarm Review

AV-Vital Vital Signs Monitor

Aveus Vital Signs Monitor provides essential measurement of patient's SpO2, NIBP and Temperature. It offers an effective solution and brings high performance vital signs monitoring at your finger tips.



Portable
Design



Infrared Ear
Thermometer

AV-Vital

Vital Signs Monitor

Wall Mount



Rolling Stand




Bed Rail
Clamp



- 5.8" high resolution color TFT-LCD screen
- Lightweight & portable design
- User-friendly interface
- Flexible configurations to meet different clinical needs
- Auto/Manual/Continuous/Average BP mode
- Real-time measurements & trend display
- Built-in rechargeable Lithium-ion battery for 8 hours of continuous work
- USB flash disk support
- Barcode scanner support
- Bi-directional communication with CMS-2000 central monitoring system
- Compatible with PatientCare Viewer PC data management software
- Wired LAN/Wi-Fi support, HL7 support via XML files or CMS-2000
- Nurse call
- Suitable for adult, pediatric and neonatal patients

Seamless
Connectivity

CMS-2000



Web Viewer

Configuration

	SpO2	NIBP	Quick Temp	Infrared Ear Thermometer
AV-Vital (SpO2 only)	✓	–	–	–
AV-Vital (NIBP only)	–	✓	–	–
AV-Vital (NIBP+Quick Temp)	–	✓	✓	–
AV-Vital (NIBP+Infrared Ear Thermometer)	–	✓	–	✓
AV-Vital (SpO2+NIBP+ Quick Temp)	✓	✓	✓	–
AV-Vital (SpO2+NIBP+ Infrared Ear Thermometer)	✓	✓	–	✓



AVEUS INSTRUMENTS Sp. z o.o.

17/19/133 Domaniewska Street,
02-672 Warsaw Poland

ul. Domaniewska 17/19 lok.
133 02-672 Warsaw Poland

info@aveusinstruments.com.pl
www.aveusinstruments.com.pl