



AVEUS



WITH
RECHARGEABLE
BATTERY

AV-WS-1

Digital Physician Scale

With Electronic Height Rod

www.aveusinstruments.com.pl

AV-WS-1

Digital Physician Scale

The Aveus **AV-WS-1** Digital Physician Scale are suitable for many medical and sports applications like gyms, doctors offices and hospitals. They feature a sturdy design and construction with a metal base and non slip weighing platform. The Aveus **AV-WS-1** is a digital scale with built in battery and color coded keypad that is simple to use.

The Aveus **AV-WS-1** scale is a Digital Physician Scale. It is an accurate electronic device with advanced design and stable performance. It is designed to measure both the weight and the height of a person. Prescribed by dieticians and physicians, the inclusive BMI (Body Mass Index) feature is easy to use for the monitoring of personal health/medical needs.

AV-WS-1 Standard Features

- Robust metal housing
- Large Anti-slip platform
- Large, clear display
- BMI-Body mass index
- Hold Function
- External calibration
- Full tare range
- Overload protection
- 2 weighing units (kg,lb)
- Rechargeable battery / mains adapter supplied as standard
- Electronic height rod as standard
- USB Interface

Specifications

- Capacity: 250kg x 0.1kg/550lbs x 0.2lb;
- Height Rod Range: 124 cm ~ 220 cm
- Graduation: 0.2cm
- LCD with backlit display
- Platform: 480 x 320 mm
- Power supply: AC Adapter
- 100-240V 50/60 Hz
- Battery: 1x lithium battery
- Continuous use cycle: 300 Hours



AVEUS INSTRUMENTS Sp. z.o.o.

17/19/133 Domaniewska Street,
02-672 Warsaw - Poland

ul. Domaniewska 17/19 lok.
133 02-672 Warszawa - Poland

info@aveusinstruments.com.pl
www.aveusinstruments.com.pl