



AVEUS

LARYNGOSCOPE MACINTOSH AV-L-5

- The laryngoscope blades are manufactured from antimagnetic 18/8- quality stainless steel type AISI 303/304 which is highly resistant to corrosion.
- AVEUS blades are maintenance free and autoclavable up to 134°C / 5 min approximately for 1,000 times.

Macintosh Blades:

- **Aveus** provide lightweight semi-rigid nylon case containing complete set of macintosh conventional 5 sizes of blades with handle.
- **Blades Size:** (0,1,2,3,4)

CE EC REP



IPXI ⚠️ ⚡️