


LARYNGOSCOPE MACINTOSH & MILLER AV-LM-5



- The laryngoscope blades are manufactured from antimagnetic 18/8- quality stainless steel type .
- AISI 303/304 which is highly resistant to corrosion. AVEUS blades are maintenance free and autoclavable up to 134°C / 5 min approximately for 1,000 times.

Blades Size

- Aveus provide 5 different sizes of blades in one packing. Macintosh Blades (2,3,4) Miller Blades (1,2)



Aveus Instruments Sp. z o.o.



17/19/133 Domaniewska Street
02-672 Warsaw – Poland

info@aveusinstruments.com.pl
www.aveusinstruments.com.pl