

# OPHTHALMOSCOPE

## AV-OP-2



- The Aveus ophthalmoscope is designed for examination of the eye.
- Made from impact resistant re-enforced brass.
- Enhanced Xenon/Halogen super bright light.
- Easy replacement of lamp.
- 5 different apertures. Large Spot, Small Spot, Semi Circle, Eye Spot, Red Free Filter



CE EC REP



IPXI   

Aveus Instruments Sp. z o.o.



17/19/133 Domaniewska Street  
02-672 Warsaw - Poland

[info@aveusinstruments.com.pl](mailto:info@aveusinstruments.com.pl)  
[www.aveusinstruments.com.pl](http://www.aveusinstruments.com.pl)