

The Aveus logo features the brand name in a bold, blue, sans-serif font. A red swoosh underline starts under the 'A' and extends to the right, ending under the 'S'.

LARYNGOSCOPE MILLER

AV-M-5

- The laryngoscope blades are manufactured from antimagnetic 18/8- quality stainless steel type AISI 303/304 which is highly resistant to corrosion. AVEUS blades are maintenance free and autoclavable up to 134°C / 5 min approximately for 1,000 times.

Miller Blades:

- **Aveus** provide lightweight semi-rigid nylon case containing complete set of MILLER conventional 5 sizes of blades with handle.
- **Blades Size:** (0,1,2,3,4)

CE EC REP



IPXI   

Aveus Instruments Sp. z o.o.
 17/19/133 Domaniewska Street
02-672 Warsaw - Poland

info@aveusinstruments.com.pl
www.aveusinstruments.com.pl