

Air-Compressing Nebulizer

AVEUS

AERONEB

Nebulizer

TECHNICAL DATA:

Power Supply: AC220V $\pm$ 22V, 50Hz $\pm$ 1Hz

Max pressure: ≥ 0.15 MPa

Max air flow: ≥ 10 L/min

Max nebulizing rate: ≥ 0.2 ml/min (by standard cup)
 ≥ 0.45 ml/min (by optional cup)

MMD: 3.9 μ m

Residual Volume: ≤ 1 ml

Operation noise: ≤ 60 dB(A)

Weight: 1.5kg

Operation time: 30minutes on, 30minutes off

G.W./N.W.: 19kg / 14.7kg

Packing: 10PCS/CTN 85*30*40cm



Aveus Instruments Sp. z o.o.



17/19/133 Domaniewska Street
02-672 Warsaw – Poland

info@aveusinstruments.com.pl
www.aveusinstruments.com.pl