

AV-VAC-C



Main Technical Specifications

- Ultimate negative pressure value: $\geq 0.075\text{MPa}$
- Regulating scope of negative pressure:
0.02MPa ultimate negative pressure value
- Air exhausting speed: $\geq 14\text{L/min}$
- Reservoir: $\geq 1000\text{ mL}$, 1pc
- Power supply: AC220V/110V and 50/60Hz
- Input power: 90VA
- Noise: $\leq 65\text{dB(A)}$
- Suction pump: piston pump
- Dimension: 335mm X 150mm X 240mm

With Battery



AV-VAC-A

AV-VAC-A/C

Portable Phlegm Suction Unit

AVEUS INSTRUMENTS Sp. z o.o.

17/19/133 Domaniewska Street,
02-672 Warsaw - Poland

ul. Domaniewska 17/19 lok.
133 02-672 Warszawa - Poland

info@aveusinstruments.com.pl
www.aveusinstruments.com.pl