

AVEUS



AV-DS-100

Integrated Diagnostic Series

A convenient, modular package that puts everything needed for fast, accurate diagnoses at arm's length. By combining a wall aneroid sphygmomanometer, direct ophthalmoscope, fiber optic otoscope, Ear thermometer and tip dispenser, Honsun puts all the basic diagnostic tool within easy reach. These systems are attractive, modular and easy to order and install. Customize your Wall Set to match your individual needs



[M1] FIBER OPTIC OTOSCOPE

- Lightweight design with high performance
- Tunable Light Source.
- 3X Magnification.
- Fiber optic Illumination
- Xenon Halogen Technology/LED Technology (optional).
- Remark: with insufflation bulb



[M3] EAR THERMOMETER

- Ear thermometer, 1 second measurement.
- Multiple use (Wide Range Measurement)
- CF Switchable
- 0 Data Memory



INSUFFLATION BULB



[M2] DIRECT OPHTHALMOSCOPE

- Lightweight design and easy to operation
- Tunable Light Source.
- 6 apertures with red-free and blue filter.
- 29 kinds of adjustable corrective range of lenses: -35D-40D.
- Xenon Halogen bulb/LED bulb (optional)



[M4] SPHYGMOMANOMETER

- Ergonomics design
- The Square shaped dial 145x145mm
- Back box is easy to hold the cuff
- High quality Sphygmomanometer



[M1.1] TIP DISPENSER (WITH DISPOSABLE TIPS)

- Tip Size: 2.75mm dia. (24pcs)
- 4.25mm dia. (24pcs)

AVEUS INSTRUMENTS Sp. z o.o.

17/19/133 Domaniewska Street,
02-672 Warsaw - Poland

ul. Domaniewska 17/19 lok.
133 02-672 Warszawa - Poland

info@aveusinstruments.com.pl
www.aveusinstruments.com.pl

