



Auto CPAP

AV-CP20

www.aveusinstruments.com.pl



Auto CPAP

AV-CP20




Integrated humidifier with optional heated air delivery system

Precise control of temperature and humidity levels




Optional remote data transfer via GPRS module



Real-time user data display

Advanced software for comprehensive sleep analysis




Enhanced trigger responsiveness for improved ventilation

Optimized airflow through reduced internal resistance



Multiple alarms to safeguard patient treatment



Optional Bluetooth oximeter for continuous SpO2 monitoring

Auto CPAP

AV-CP20

Technical Specifications

Product model	Auto CPAP (AV-CP20)
Work mode	CPAP AUTO
Pressure Range	4-20 cmH ₂ O, automatic adjustment
Expiratory pressure	-
Inspiratory pressure	-
Tubing heating function	Yes
Data cloud platform	Optional
Mobile network	Optional
Pressure rise/fall time	-
Pressure difference	-
I/E sensitivity	-
Ramp time	0-60min
Noise	≤ 22dBA @10cmH ₂ O
Data storage	TF card, >5 years
Screen (inch)	3.5inch color display
Humidity level	Heating and humidification, level 1-5 adjustable
Heating tubing temperature	10°C-30°C adjustable, step size 1°C
Preheat	ON/OFF
Flight mode	ON/OFF
Bluetooth oximeter	Optional
Alarm	Mask off, suffocation, pressure abnormality, tubing abnormality
Power	Input: 100-240V~, 50/60Hz, 1.5A Max. Output: 24V - 3A
Dimensions(cm)	26.68Lx17.25Wx11.85H(machine)
Carton dimensions(cm)	41L x 29W x 17.5H (1pc/carton)
	41L x 32.5W x 30H (2pcs/carton)
	48L x 41W x 30H (3pcs/carton)
Weight(with humidifier)	4.07Kg (1pc/carton)
	7.74Kg (2pcs/carton)
	11.57Kg (3pcs/carton)

AVEUS INSTRUMENTS Sp. z o.o.

17/19/133 Domaniewska Street,
02-672 Warsaw - Poland

ul. Domaniewska 17/19 lok.
133 02-672 Warszawa - Poland

info@aveusinstruments.com.pl
www.aveusinstruments.com.pl

