



AV-ML-1

Oxygen Minilator

AV-ML-1

Oxygen Minilator

Detailed Specifications

Model	AV-ML-1
Supply Pressure	40-90 psig
Operational Temperature	-30° F to 125° F
Storage Temperature	-40° F to 160° F
Flow Rate	4, 6, 10 OR 15 lpm
Materials	Manifold, Fittings, Barbs, Orifices
Flow Rate	4, 6, 10 OR 15 lpm

Product Information :

Oxygen Minilator is a versatile constant flow device designed to provide oxygen for up to five patients simultaneously.

Featuring durable anodized aluminum construction and 50 psig outlets, it supports various medical devices and flow rates (4, 6, 10, or 15 lpm). Ideal for emergency and mass casualty situations.



Expandable for
Emergency Situation




Compact and
Lightweight Design



Durable Anodized
Aluminium



2M Reusable Hose with
DISS Adapter at both sides



AVEUS INSTRUMENTS Sp. z o.o.

17/19/133 Domaniewska Street,
02-672 Warsaw - Poland

ul. Domaniewska 17/19 lok.
133 02-672 Warszawa - Poland

info@aveusinstruments.com.pl
www.aveusinstruments.com.pl