



Portable X-Ray

AV-XR-1

AV-XR-1 is one Multi-functional Portable X-ray System, featured by small size, light weight, high power density by 5-KW High Voltage Generator, fit for all type portable X-ray unit applic


Features



Small & Portable, Easy to carry



Carry Wireless Exposure



Cable Exposure



Easy Operation, Fast Imaging


Support Bluetooth Connection to Laptop



Portable X-Ray

AV-XR-1

Application Way




Portable X-Ray

AV-XR-1

Technical Specifications

Model	Portable X-Ray (AV-XR-1)
Power Condition	
Power Input	AC 210~230V/AC 100V~120V
Power Rate	5KW
High Voltage Generator	
kV Range	40~125 kV
mAs Range	0.1~100mAs
mA Range	10~100mA
Communication	RS232/Bluetooth
Exposure	Wired/Wireless
X-ray Tube	
Dual Focus	0.6/1.8mm
Target Angle	15°
Physical Specification	
Size	505 x 348 x 230mm
Weight	16kg



AVEUS INSTRUMENTS Sp. z o.o.

17/19/133 Domaniewska Street,
02-672 Warsaw - Poland

ul. Domaniewska 17/19 lok.
133 02-672 Warszawa - Poland

info@aveusinstruments.com.pl
www.aveusinstruments.com.pl

