

POWER WHEELCHAIR

ENVOL - AVEUS RYDER



One-touch fold and unfold design, Detachable battery design, Lightweight & Portable, Accurate and sensitive controller, PU solid tire, Electromagnetic braking system.

Fold Size	590X275X760mm
Unfold Size	950X590X940mm
Material	Carbon Fibre
Max Weight Capacity	140kg
Speed	0-6Km/h
Operating Range	20Km
Max. Inclination	8°
Turning Radius	80cm
Motor	24v 180w
Batteries	24v 6AHX2 Ternary Lithium Battery
Battery Charge	24v 2A
Front Wheel Size	7"
Rear Wheel Size	8"
Seat Width	400mm
Weight	17Kg



Aveus Instruments Sp. z o.o.



17/19/133 Domaniewska Street
02-672 Warsaw – Poland

info@aveusinstruments.com.pl
www.aveusinstruments.com.pl