



Multi-Functional Electrosurgical Unit

AV-ESU-400

INTRODUCTION

High-frequency Electrosurgical Unit is an advanced electro surgery equipment for incision and coagulation of the biological tissue through the high frequency current.

It is an excellent substitute for mechanical scalpels due to its simple & safe operation, fast tissue incision and blood coagulation. It is wide use in General surgery, Gynecology, Urology, Oncology, Orthopedics etc.



Multiple Functions

Both monopolar and bipolar work modes.



Easy-to-Use

Simple buttons to convert the work modes and change the power output.



Clear Display

3 digital electronic LED for the power output and error codes.



High Safety

REM (Return Electrode Monitoring) system to avoid burnt or damage to the surgeon and patient.



Wide Application

It can be connected to endoscopes for more clinical applications.

Features:

- Easy to use buttons for simple operation;
- High power output 400W for more operation selection;
- Up to 7 work modes: Pure/Blend1/Blend2, Standard/Spray,Bipo1/Bipo2
- Increment of 5W micro power control system to ensure optimal effect;
- Minimal sticking & carbonization, reduce the surgery damage;
- Auto detection of impedance to protect the patient
- Support underwater cutting (TURP etc).
- OEM projects welcome



Specification

Function		Work Mode	Frequency	Rated Power	Gain Interval	Form of HF Current
Monopolar	Cut	Pure Cut Blend 1 Blend 2	500KHz 500KHz 500KHz	0-400W/300Ω 0-250W/300Ω 0-200W/300Ω	5 Watt	Sine Wave Modulate Sine Wave Modulate Sine Wave
	Coag	Standard Spray	500KHz 500KHz	0-100W/300Ω 0-70W/300Ω	5 Watt	Modulate Sine Wave Modulate Sine Wave
Bipolar		Bipo Cut Bipo Coag	500KHz 500KHz	0-100W/50Ω 0-70W/50Ω	5 Watt	Sine Wave Modulate Sine Wave

AV-ESU-400

Power Supply:

- 110V, 50/60Hz or 220V,50/60Hz

Dimension:

- Product: 421(L) x 357(W) x 162(H) mm
- Package: 520(L) x 445(W) x 305(H) mm

Weight:

- Net Weight: 8.7KG
- Gross Weight: 14.5KG

AVEUS INSTRUMENTS Sp. z o.o.

17/19/133 Domaniewska Street,
02-672 Warsaw - Poland

ul. Domaniewska 17/19 lok.
133 02-672 Warszawa - Poland

info@aveusinstruments.com.pl
www.aveusinstruments.com.pl

Standard Accessory



Standard Electrodes



Push Button Electrosurgical Pencil



Foot Switch Electrosurgical Pencil



Foot Switch



Disposable Bipolar Neutral Pad



Patient Plate Cable

Optional Accessory



Bipolar Forceps



Silicon Patient Plate



Metal Patient Plate

AVEUS®