

AURA CROSS

Emergency Stretcher

The AURA CROSS by AVEUS is a compact and lightweight self-loading stretcher designed for patient transport. Ideal for rescue operations in both road and medical vehicles, its legs fold independently using dedicated controls and automatically extend during unloading. This model offers a higher loading capacity.



SPECIFICATIONS

- Aluminium board easy to clean
- Adjustable backrest
- 2 fixed wheels and 2 pivoting with large diameter to make movement easier
- Brakes on rear wheels
- Foldable side bar Safe bar
- Adjustable heights system

TECHNICAL DATA

Length	1900±10mm
Width	65±10mm
Weight	46±1 kg
Loading Capacity	300kg
Wheels	φ 200mm
Height of the loading carriage wheel	63.5 & 66.5cm
Side Bars Length	650mm
Side Bars Height	220mm

AVEUS INSTRUMENTS Sp. z.o.o.

17/19/133 Domaniewska Street,
02-672 Warsaw - Poland

info@aveusinstruments.com.pl
www.aveusinstruments.com.pl