

AURA TEC

Emergency Stretcher

The AURA TEC by AVEUS is a durable, adjustable medical transport device with an aluminum board, safety bars, and maneuverable wheels. It features an adjustable backrest, height system, and a 350kg loading capacity, ensuring safe and efficient patient handling.



SPECIFICATIONS

- Aluminium board easy to clean
- Adjustable backrest
- 2 fixed wheels and 2 pivoting with large diameter to make movement easier
- Brakes on rear wheels
- Foldable side bar Safe bar
- Adjustable heights system

TECHNICAL DATA

Length	2000±10mm
Width	65±10mm
Weight	45±1 kg
Loading Capacity	350kg
Wheels	φ 200mm
Height of the loading carriage wheel	61, 65 & 69cm
Side Bars Length	840mm
Side Bars Height	220mm

AVEUS INSTRUMENTS Sp. z.o.o.

17/19/133 Domaniewska Street,
02-672 Warsaw - Poland

info@aveusinstruments.com.pl
www.aveusinstruments.com.pl